

THE LOOMING PROBLEM OF **RESISTANCE**



Faecal Worm Egg Counts using **FLOTAC Technique**

A new, more sensitive* counting method also giving early warning of resistance

Only commercially available
in the UK from
Langford Veterinary Services

Special Price
Only £10
per sample

This highly sensitive test, using the **FLOTAC technique** (*a new improved counting method developed in Italy that is accurate to 1 egg per gram (epg) compared to the McMaster method, 25-50 epg) will not only give you very accurate

egg counts but also alert you to an early warning of resistance.

If you have clients whose animals would benefit from this test, please look at our website or contact us directly for further information.



Langford

VETERINARY SERVICES

- t** +44 (0)117 928 9412
- f** +44 (0)117 928 9613
- e** lvslab@langfordvets.co.uk
- w** www.langfordvets.co.uk

a **Langford Veterinary Services**, Diagnostic Laboratories
Churchill Building, Langford House, Langford, Bristol BS40 5DU

Langford Veterinary Services is a fully owned subsidiary of the University of Bristol